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Water based  
Surface Disinfection

**FIELD OF APPLICATION**

VERIFORTE® is a rapid and effective water based disinfectant for use on surfaces. VERIFORTE® is pH neutral, alcohol-free, aldehyde-free, fragrance-free and non flammable.

**INSTRUCTIONS FOR USE**

VERIFORTE® is a ready to use, fluid, clear solution. Surfaces and objects are simultaneously und thoroughly moistened with 50ml/m<sup>2</sup> using where appropriate a disposable and leaving to dry. Rinsing with water is not necessary. The product evaporates leaving no residue.

**HIGHLY EFFECTIVE**

VERIFORTE® is tested according to European Standards EN 13697, EN 13704 and EN 14476 and is listed by DGHM/VAH. VERIFORTE® can be used both as a wipe and spray disinfectant. The stated exposure times below, demonstrate the full effectiveness of Veriforte against bacteria, viruses, spores and yeast, even when, mechanical help is not employed.

**EXPOSURE TIMES** undiluted without mechanical help:

<b>Bacterial Effect</b>	surface precleaned 5min*	
	MRSA.....5 min	E.coli.....1 min
	E.hirae.....1 min	P.aeruginosa.....1 min
<b>Yeasticidal Effect</b>	precleaned surface	
	C. albicans.....5 min	
	*shortest possible exposure time according to VAH – Applied Hygiene Network	
<b>Sporocidal Effect</b>	precleaned surface	
	Bac. Subtilis.....5 min	
<b>Virucidal Effect</b>	Murine Norovirus	low burden 1 min / high burden 1 min
	Adenovirus	low burden 1 min / high burden 1 min
	Poliovirus	low burden 1 min / high burden 30 min
<b>Virucidal against</b>	enveloped viruses.....5 min according to the DVV and RKI guidelines	
(incl. HIV, HBV und HCV)		

**SPECIFICATIONS** Color/Form: transparent/liquid pH value: 6.5---7.5

**COMPOSITION** >99.5% Water; <0.4% NaCl (Sodium Chloride); <0.05% resp. 0.5g/L NaOCl (Sodium Hypochlorite) / HOCl (Hypochlorous acid)

**REGISTRATION**  
 N – 54547 for Product type 2 (disinfectants)  
 N – 54587 for Product type 4 (Food and feed sector)  
 N – 54588 for Product type 5 (drinking water)

**HAZARDS IDENTIFICATION**

The product is non-hazardous. Normal precautions should however be observed when handling chemicals: Note the lot number and expiration date printed on the packaging.

**SAFETY**

In normal use, VERIFORTE® is not dangerous for people. Avoid eye contact. Irritation may occur. After eye contact, rinse immediately with clear water for several minutes (remove contact lenses). If you experience lengthy complaints related to the product, please consult a medical specialist.

**STORAGE**

Store VERIFORTE® in a light-protected and cool place.  
 Keep container closed tightly. Note the lot number and expiration date printed on the packaging.

**DISPOSAL**

With a sufficient dilution Veriforte can be disposed via a standard wastewater sewage system. Please observe local and country---specific waste and wastewater regulations. After emptying and thorough rinsing with water, the plastic container can be recycled or can be disposed of in an appropriate way.

**SYMBOLS**

 Temperature limits

 Protect from sunlight

 Period after opening in months

 Read product information before use

**Use biocides safely. Always read the label and product information before use.**

Distribution: CSE Clean Solution GmbH | Salzachtalstraße 86 | A-5400 Hallein | T: +43 (0) 6245 75 642-0 | office@cse-cleansolution.com

The information in this product data sheet corresponds to our current level of knowledge. They are the result of our previous experience and scientific reviews. However, we can not individually examine and influence different circumstances in practical application. The quality of the work therefore depends on the professional processing by the user. The user is obligated to always check the compatibility of the material to be treated with the product in an inconspicuous place or to obtain technical advice. Liability for incomplete or incorrect information in our data sheet exists only in case of gross negligence or intent. Any claims arising from the Product Liability Act remain unaffected. With this product data sheet, all previous product information will lose their validity. Document Name: DI-7-205-2, February 2020